FIBRI ONE

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Telefone: +55 11 40231749 | e-mail: contact@biocontinental.com

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- 1- Hirsh J., Deykin D., Poller L.: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2- International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
- 3- National Committee For Clinical Laboratory Standards: Colection, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4 - Douketis Jd, Lane A, Milne J, Ginsberg Js. Trombosis Research, 1998; 92: 11-17.

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FIBRI ONE

CLINICAL SIGNIFICANCE

Fibrinogen is a glycoprotein present in plasma and platelet granules, this being the coagulation factor (Factor I) that is found in the highest plasma concentration (150 - 350 mg/dL). Fibrinogen originates in the liver and is converted to fibrin during coagulation. Because fibrin is necessary for the formation of clots, fibrinogen deficiency can produce mild to severe bleeding disorders. Also known as acute phase protein, fibrinogen can occur at high levels in patients with acute infection, collagen disorders, neoplasms, nephrosis or hepatitis. Elevated levels can also occur during the third trimester of pregnancy and in women taking oral contraceptives. Depressed levels of fibrinogen may indicate congenital afibrinogenemia, hypofibrinogenemia, DIC, fibrinolysis, severe liver disease, metastatic cancer or bone marrow lesions. Trauma or obstetric complications can also cause low levels. Fibrinogen levels below 100 mg/dL prevent an interpretation of coagulation tests that have a fibrin clot as an end point.

PURPOSE AND BASIS OF THE METHOD

Method for quantitative determination of fibrinogen in citrated plasma. Only for *in vitro* diagnostic use. It consists of the quantitative determination of fibrinogen in citrated plasma, after the addition of a standardized amount of thrombin in a sample of diluted citrated plasma and measurement of the clotting time. When the thrombin concentration is high, the clotting time of the diluted citrated plasma is inversely proportional to the concentration of fibrinogen. To obtain the fibrinogen concentration, the analyzed plasma clotting time is compared with the plasma clotting times with known fibrinogen concentration.

REAGENTS PROVIDED

Number 1 - Buffer - Store between 2 and 10°C. Imidazole > 10mmol / L and preservative <1%.

Number 2 - Thrombin - Store between 2 and 10°C. Thrombin 100 NIH Units / mL, Albumin > 1%, buffer > 10 mmol / and preservative > 1%.

Number 3 - Reference plasma - Store between 2 and 10°C. Buffer > 25 mmol / L and lyophilized plasma. Potentially infectious.

Note: See concentration on the bottle label.

MATERIALS REQUIRED AND NOT PROVIDED

Coagulometer, 37°C water bath, clock or timer, pipettes, test tubes, and control plasmas. These items are found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 10°C. Transport, at temperatures up to 30°C, should not exceed 5 days. Keep away from light and avoid humidity.

PRECAUTIONS

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the proposed methodology to obtain exact results.

3- The water level in the water bath must be higher than the level of the reagents in the test tubes.

4- Avoid heating the reagent that will not be used. Reagent that is already heated should not be returned to the original bottle.

5- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the Product Safety Data Sheet.

7- Do not use the product in case of damage to the packaging.

8- It is essential that the equipment used is properly calibrated and subjected to periodic maintenance.

SAMPLES

Material: Plasma

Anticoagulant: Use Sodium Citrate 3.8% (0.130 M) or 3.2% (0.109 M) as anticoagulant, in the proportion of 9 parts of Blood to 1 of Anticoagulant (as an example: 4.5 mL of Blood + 0.5 mL of Anticoagulant), or commercial Citrate kits, according to the manufacturer's Instructions for Use, or vacuum tubes with citrate, available on the market.

Collection: Blood must be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, bubble formation and tissue fluid aspiration. When evaluating only the coagulation test, always collect 2 tubes of citrate and discard the first.

Preparation: The blood should be mixed with the anticoagulant immediately after collection.

Centrifuge immediately at 3000 rpm for 15 minutes.

Remove the plasma without pipetting red cells and the yellow layer.

The samples must be tested in less than 3 hours. If the test cannot be performed during this period, the plasma should be frozen for a maximum of 1 week at -20°C.

POST OPEN STABILITY

The stability of the FIBRI ONE kit after opening is 14 days.

PROCEDURE

Plasma Preparation Reference:

Carefully open the bottle and add 1 mL of deionized water. Homogenize. Leave the bottle to stand for 30 minutes. Every 10 minutes carefully shake the bottle with circular movements to prevent the material from adhering to the wall of the bottle. Avoid foaming.

Preparation of dilution for calibration curve, samples and control:

1- Prepare the dilutions for mounting the calibration curve according to table 1 and perform the technique from item 3, considering each dilution as a sample.

2- Dilute sample and 1:10 control (0.1 mL of sample + 0.9 mL of R1 Buffer). Perform the technique starting from item 3

TEST

3- Incubate 0.2 mL of the diluted sample in a water bath at 37°C for 4 to 6 minutes.

4- Add 0.1 mL of the unheated R2 (Thrombin), simultaneously start the timer and record the clot formation times.

CALIBRATION CURVE

Duplicate clotting times for each dilution R3 (Plasma Reference) and obtain the average time. Build a graph on a log-log paper applying the concentration values on the X axis (abscissa) and the clotting times on the Y axis (ordered).

After knowing the patient's and control's clotting time, place the result on the calibration curve and read the fibrinogen value for the time found.

Table 1 - Calibration curve

R3 - Plasma Reference		R1 Buffer (mL)	Dilution Factor	Fibrinogen = Concentration R3 (mg/dL) x Dilution Factor	Time (sec)
Dilution	Volume (mL)				
1:5	0.1	0.4	10/5 = 2.0	- X 2.0	-
1:10	0.1	0.9	10/10 = 1.0	- X 1.0	-
1:15	0.1	1.4	10/15 = 0.67	- X 0.67	-
1:20	0.1	1.9	10/20 = 0.5	- X 0.5	-
1:30	0.1	2.9	10/30 = 0.33	- X 0.33	-

The R3 - Plasma Reference diluted 1:10 represents 100% of the value indicated on the bottle label, that is, the value indicated on the label corresponds to the concentration of the plasma already diluted 1:10.

RESULTS

If the clotting time of the sample is less than the clotting time of the 1:5 dilution of the curve, repeat the test by diluting the sample 1:30, read the result on the calibration curve and multiply it by 3.

If the clotting time of the sample is greater than the clotting time of the 1:30 dilution of the curve, repeat the test by diluting the sample 1:3, read the result on the calibration curve and multiply it by 0.3.

If the 1:3 dilution clotting time is greater than the longer clotting time of the calibration curve, the result should be reported as less than the lowest concentration of the calibration curve.

TRACEABILITY

The Plasma Reference of the FIBRI ONE kit is traceable to the reference material NIBSC 98/612 (International Standard for Fibrinogen Plasma - WHO International Standard).

REFERENCE VALUES

The range should be used as a guide only. It is recommended that each laboratory establishes, in the population served, its own reference intervals, using a pool of normal plasmas or a reference plasma.

Fibrinogen: 150 to 350 mg / dL

Conversion: Conventional units (mg / dL) x 0.01 = SI units (g / L)

INTERFERENT

Icemic, lipemic or hemolyzed samples.

Use of plasmas containing Heparin, EDTA or Sodium Oxalate.

Elevated levels of fibrinogen or fibrin degradation products can prolong clotting times especially when fibrinogen levels are less than 150 mg/dL.

PRODUCT PERFORMANCE

COMPARISON OF METHODS

The FIBRI ONE Kit was compared with another commercially available kit for determining fibrinogen. Seven analyzes were performed and the results were evaluated. The linear regression coefficient: $y = 1.0082x - 3.4767$ and the correlation coefficient = 0.9992. With these results, it can be concluded that the kit has good methodological specificity.

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Concentration Average (mg/dL)	131.2	231.75	311.85
Standard Deviation (sec)	1.11	1.71	1.66
Coefficient of Variation (%)	0.84	0.74	0.53

REPRODUCIBILITY

Reproducibility was calculated from 20 successive determinations over 3 consecutive days, using 3 samples with different concentrations, obtaining the following results::

	Sample 1	Sample 2	Sample 3
Concentration Average (mg/dL)	131.9	232.317	312.867
Standard Deviation (sec)	1.80	2.22	2.49
Coefficient of Variation (%)	1.37	0.96	0.79

PRESENTATION

Presentación	Reactivo 1	Reactivo 2	Reactivo 3
1	1 x 60 mL	5 x 1 mL	1 mL
2	2 x 60 mL	10 x 1 mL	1 mL
3	4 x 60 mL	20 x 1 mL	1 mL
4	60 mL	5 mL	1 mL
5	2 x 60 mL	1 x 10 mL	1 mL
6	2 x 60 mL	2 x 5 mL	1 mL
7	4 x 60 mL	2 x 10 mL	1 mL
8	4 x 60 mL	4 x 5 mL	1 mL

CONSUMER INFORMATION:

Biocontinental guarantees the quality of the product, provided that all conditions are respected, and the storage and transport care, specified in the instruction for use, until the expiration date indicated on the packaging.

If further information or guidance is required, the customer should contact Biocontinental Biological Products.

Technical Assistance and request for Instruction of Use in printed format:

Phone: +55 11 40231749 | E-mail: contact@biocontinental.com

BIBLIOGRAPHIC REFERENCES

1 - Hirsh J, Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.

2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88: 779 1985.

3 - National Committee For Clinical Laboratory Standards: Colleccon, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.

4 - Douketis Jd, Lane A, Milne J, Ginsberg Js. Trombosis Research, 1998; 92: 11-17.

Manufacturer | Regularized by: Continental Produtos Biológicos Ltda
 Address: Rua Santana, 305, 13300-220 Itu- SP
 Site: www.biocontinental.com
 E-mail: contact@biocontinental.com
 CNPJ: 61.058.400/0001-42
 Brazilian Industry

ANVISA: 80132390001

Attention user: check if the version obtained from the Instruction for Use not printed corresponds to the product purchased.

Version: JUN/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED