## THROMBO ST

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- 1- Hirsh J, Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2- International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
- 3- National Committee For Clinical Laboratory Standards: Collection, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4- Douketis Jd, Lane A, Milne J, Ginsberg Js. Trombosis Research, 1998; 92: 11-17.

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Versão: MAR/2025

## UNIVERSAL SYMBOLOGY

|                                                                                       |                                       |                                                                                       |                                        |
|---------------------------------------------------------------------------------------|---------------------------------------|---------------------------------------------------------------------------------------|----------------------------------------|
|    | CATALOG NUMBER                        |    | MANUFACTURED BY                        |
|    | BATCH CODE                            |    | CONTROL                                |
|    | DATE OF MANUFACTURE                   |    | POSITIVE CONTROL                       |
|    | USED BY<br>(last day of month)        |    | NEGATIVE CONTROL                       |
|    | TEMPERATURE LIMITATION<br>(store at)  |    | BIOLOGICAL RISK                        |
|   | CONTAINS SUFFICIENT<br>FOR <N> TESTS  |   | INFLAMMABLE                            |
|  | CONSULT INSTRUCTIONS<br>FOR USE       |  | CORROSIVE                              |
|  | IN VITRO DIAGNOSTIC DEVICE            |  | POISON                                 |
|  | EUROPEAN AUTHORIZED<br>REPRESENTATIVE |  | CE MARK                                |
|  | KEEP AWAY<br>FROM SUNLIGHT            |  | DO NOT USE IF<br>PACKAGE IS<br>DAMAGED |

## THROMBO ST

### CLINICAL SIGNIFICANCE

Prothrombin time or Quick time is the most clinically important screening test in the evaluation of disorders of the extrinsic coagulation pathway. The test measures the Prothrombin Time, this increased time is due to congenital or acquired deficiencies of factors I, II, V, VII and X. Being thus, the test of choice for assessing the reduction of the factors vitamin K dependent, and indicated in the control patients using oral anticoagulants, hemorrhagic disease of the newborn, hepatic failure, intravascular coagulation, disordered eating and absorption of vitamin K. Evaluates both the extrinsic pathway as the common pathway of coagulation.

### PURPOSE AND BASIS OF THE METHOD

Method for determination of manual or automated of Prothrombin Time in citrated plasma. For *in vitro* diagnostic use only. Consists of measuring the Clotting Time of plasma after addition a source of Thromboplastin (factor III, tissue factor) and Calcium. The factor VII is then activated to form a complex (complex FT-FVIIa) that will activate factors IX and factor Xa XO with Phospholipids factor tissue, Calcium and factor Va form the Prothrombin activator complex which transforms Prothrombin (factor II) into Thrombin (factor IIa). Thrombin acts on Fibrinogen, transforming it into Fibrin. This test assesses the factors of Prothrombin Complex (factors II, V, VII and X).

### REAGENTS PROVIDED

**Number 1 - Clot Former** - Store between 2 and 8°C. Contains: Rabbit brain extract and preservative.

### Presentation ready for use.

### MATERIALS REQUIRED AND NOT PROVIDED

Log-log paper, pipette and tips, tubes, coagulometer or stopwatch and plasmas controls. Can be found in stores specialized in articles for Laboratories of Clinical Analysis.

### STORAGE, STABILITY AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 10°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep out of the light. Do not freeze.

### PRECAUTIONS

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The water level in the water bath must be higher than the level of the reagents in the test tubes.
- 4- Avoid heating the reagent that will not be used. Reagent that is already heated should not be returned to the original bottle.
- 5- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 6- To obtain information related to biosafety or in case of accidents with the product, consult the Product Safety Data Sheet.
- 7- Do not use the product in case of damage to the packaging.
- 8- It is essential that the equipment used is properly calibrated and subjected to periodic maintenance.

### SAMPLES

**Material:** Plasma

**Anticoagulant:** Use Sodium Citrate 3.8% (0.130 M) or 3.2% (0.109 M) as anticoagulant, in the proportion of 9 parts of Blood to 1 of Anticoagulant (as an example: 4.5 mL of Blood + 0.5 mL of Anticoagulant), or commercial Citrate kits, according to the manufacturer's Instructions for Use, or vacuum tubes with citrate, available on the market.

**Collection:** Blood must be obtained through veinpuncture, avoiding hemolysis, prolonged tourniquet, bubble formation and tissue liquid aspiration. When evaluating only the coagulation test, always collect 2 tubes of citrate and discard the first one

**Preparation:** Blood should be mixed with anticoagulant right after collection. Centrifuge immediately at 3000 rpm, during 15 minutes.

Remove plasma without pipetting red cells and the yellow layer.

Samples must be tested in less than 3 hours. If tests can not be made during this period plasma should be frozen for at maximum 1 week at -20°C.

### POST OPEN STABILITY

The calibration stability of the THROMBO ST kit after opening is 14 days. This stability may vary according to the test, equipment and environment conditions. Therefore, it is suggested to monitor the product's performance using control plasmas.

### PROCEDURE

- 1- Preheat the Reagent N° 1 (Clot Former) at 37°C for 4 to 5 minutes (Heat only what is necessary for the test).
- 2- Homogenize the reagent before use.
- 3- In a clean tube, pipet 100 µL of the plasma to be measured (Controls or Samples) and incubate at 37 °C for 2 to 3 minutes (in a water bath or thermal blocks).

4- Add 200 µL of reagent N° 1 (Clot Former), previously heat, to the tube containing the plasma and trigger the timer.

5- Leave the tube on the heating wand slowly.

6- Remove the tube before the expected coagulation time (approximately 9 seconds), observe the formation of the clot. Immediately stop the timer by observing the formation of the clot and recording the time.

### CALCULATIONS

The results of Prothrombin Time may be expressed as:

**Quick Time:** In seconds.

**Percentage of activity:** Related to a normal plasma (considered 100% activity), and for this you must draw a curve activity, using a normal plasma pool and freshly harvested and diluted with saline (0.9%).

### CURVE

| Dilution                  | 1:1 | 1:2 | 1:3 | 1:4 | 1:8  | 1:10 |
|---------------------------|-----|-----|-----|-----|------|------|
| Pool (mL)                 | 0.2 | 0.2 | 0.2 | 0.2 | 0.2  | 0.2  |
| Physiologic Solution (mL) | —   | 0.2 | 0.4 | 0.6 | 1.4  | 1.8  |
| Activity (%)              | 100 | 50  | 33  | 25  | 12.5 | 10   |

In a log-log paper, plot the graph by putting the seconds value on the "y" axis and the percentage value on the "x" axis. This procedure should be repeated on every batch.

### INR (International Normalized Ratio)

To obtain INR must be related to the prothrombin time ratio the patient by the time of normal plasma with the ISI (Sensitivity Index International), as shown below. **The value of ISI is sent with the kit.**

### Calculation

INR = (PT patient / PT normal) ISI

Patient PT: 13.5 sec

Normal PT (reference): 12 sec

ISI: 1.20

INR = (13.5 / 12) 1.20 = 1.15

### TRACEABILITY

The THROMBO ST kit is traceable to the reference material NIBSC RBT/05 (International Standard for Thromboplastin - WHO International Standard).

### REFERENCE VALUES

**Prothrombin time (quick)** - 10 to 14 seconds

**Activity (%)** - 70 to 100%

**INR** - 1.0 to 1.08 in healthy people

- 2.0 to 3.5 in patients taking oral anticoagulants

Each laboratory should establish its reference values, and the results vary with respect to the reagent used. The results should be expressed in INR to normalize the variation of results between laboratories.

ISI values are determined by comparison with a reference material for thromboplastin primary. As low as the ISI is, the more sensitive is the reagent.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

### INTERFERENCES

Excessive hemolysis and lipemia. Use of plasma with Heparin and EDTA.

The PT can be increased by substances such as Anticoagulants Heparin and Warfarin, Oral Contraceptives, Ethanol, Tetracycline, Corticosteroids and EDTA.

The reduction of PT is observed in patients taking: Vitamin K, Anti-Histamines, Caffeine and Phenobarbital.

### PRODUCT PERFORMANCE

### COMPARISON OF METHODS

The THROMBO ST kit was compared with another method for evaluating commercially available Prothrombin Time. 25 analyzes were performed and the results were evaluated. The linear equation obtained was  $y = 1.0276x - 0.481$  and the correlation coefficient equal to 0.977. With these results, it can be concluded that the kit has good methodological specificity.

### REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 |
|------------------------------|----------|----------|
| Time (sec)                   | 14.12    | 11.87    |
| Standard Deviation (sec)     | 0.29     | 0.18     |
| Coefficient of Variation (%) | 2.08     | 1.54     |

### REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 1 sample, obtaining the following results:

|                              | Sample 1 |
|------------------------------|----------|
| Time (sec)                   | 12.477   |
| Standard Deviation (sec)     | 0.187    |
| Coefficient of Variation (%) | 1.50     |

### SENSITIVITY

The analytical sensitivity of the THROMBO ST kit is given by the value of Sensitivity Index International (ISI) and ranges from 1.2 to 1.5.

### SPECIFICITY

The specificity of the kit is provided by the addition of an exclusive source thromboplastin, specific activator of factor VII, which eliminates the possibility of direct activation of the cascade of other factors coagulation.

### PRESENTATION

| Presentation | Reagent 1 |
|--------------|-----------|
| 1            | 10 x 5 mL |
| 2            | 20 x 2 mL |
| 3            | 10 x 2 mL |
| 4            | 20 x 5 mL |

### CONSUMER INFORMATION:

Biocontinental guarantees the quality of the product, provided that all conditions are respected, and the storage and transport care, specified in the instruction for use, until the expiration date indicated on the packaging.

If further information or guidance is required, the customer should contact Biocontinental Biological Products.

### Technical Assistance and request for Instruction of Use in printed format:

Phone: +55 11 40231749 | E-mail: [contact@biocontinental.com](mailto:contact@biocontinental.com)

### BIBLIOGRAPHIC REFERENCES

- 1- Hirsh J, Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. 89/12 Chest (Suppl) 115-155, 1986.
- 2- International Committee For Standardization In Haematology International Committee On Thrombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
- 3- National Committee For Clinical Laboratory Standards: Collection, Transport And Processing Of Blood Coagulation Specimens For Testing And Performance Of Coagulation Assays, 1991. NCCLS.
- 4- Douketis JD, Lane A, Milne J, Ginsberg JS. Thrombosis Research, 1998, 92: 11-17.

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Brazilian Industry

ANVISA: 80132390004

Attention user: check if the version obtained from the Instruction for Use not printed corresponds to the product purchased.

Version: MAR/2025

## UNIVERSAL SYMBOLOGY

|                                                                                       |                                       |                                                                                       |                                        |
|---------------------------------------------------------------------------------------|---------------------------------------|---------------------------------------------------------------------------------------|----------------------------------------|
|    | CATALOG NUMBER                        |    | MANUFACTURED BY                        |
|    | BATCH CODE                            |    | CONTROL                                |
|    | DATE OF MANUFACTURE                   |    | POSITIVE CONTROL                       |
|    | USED BY<br>(last day of month)        |    | NEGATIVE CONTROL                       |
|    | TEMPERATURE LIMITATION<br>(store at)  |    | BIOLOGICAL RISK                        |
|   | CONTAINS SUFFICIENT<br>FOR <N> TESTS  |   | INFLAMMABLE                            |
|  | CONSULT INSTRUCTIONS<br>FOR USE       |  | CORROSIVE                              |
|  | IN VITRO DIAGNOSTIC DEVICE            |  | POISON                                 |
|  | EUROPEAN AUTHORIZED<br>REPRESENTATIVE |  | CE MARK                                |
|  | KEEP AWAY<br>FROM SUNLIGHT            |  | DO NOT USE IF<br>PACKAGE IS<br>DAMAGED |